


## T-Nanocrystalline Sheets (Die-cutting)

MS700-T57. 6X56. 6X0. 11A



### CORE SETS

DESCRIPTION	DIMENSIONS							$\mu'$ @128kHz	WEIGHT (g)
	A (mm)	B (mm)	T (mm)	PET ( $\mu$ m)	PS (inner) ( $\mu$ m)	Nano- (Layers)	PSA (outer) ( $\mu$ m)		
MS700-T57. 6X56. 6X0. 11A	$57.60 \pm 0.10$	$56.60 \pm 0.10$	$0.11 \pm 0.01$	$10 \pm 2$	$5 \pm 1$	4	$10 \pm 2$	$700 \pm 15\%$	2.6
MS1000-T57. 6X56. 6X0. 08A	$29.95 \pm 0.10$	$24.95 \pm 0.10$	$0.08 \pm 0.010$	$10 \pm 2$	$5 \pm 1$	3	$10 \pm 2$	$1000 \pm 20\%$	2.1
MS1400-T57. 6X56. 6X0. 06A	$29.95 \pm 0.10$	$24.95 \pm 0.10$	$0.06 \pm 0.010$	$10 \pm 2$	$5 \pm 1$	2	$10 \pm 2$	$1400 \pm 20\%$	1.6
MS2000-T57. 6X56. 6X0. 04A	$29.95 \pm 0.10$	$24.95 \pm 0.10$	$0.04 \pm 0.010$	$10 \pm 2$	/	1	$10 \pm 2$	$2000 \pm 20\%$	1.1

#### Notes :

1. The type of tape can be selected according to customer's requirements;
2. PET : Single-adhesive Tape ,PSA: Adhesive Tape;
3. The thickness of the PET&PSA is only for reference.